

**Technical Service Contact Information**

**Phone:** +98 21 86703137

**Email:** tpr.assist@gmail.com

**Hours:** S-T 8:00 AM to 3:30 PM



**TEB PAZHOHAN RAZI**

**(TPR)**

**Superoxide Dismutase Activity  
Assay Kit (48/96 Tests)**

**Storage and Stability**

This kit will perform as specified if stored at 4°C

Use before the **expiration date**

## Technical Service Contact Information

Phone: +98 21 86703137

Email: tpr.assist@gmail.com

Hours: S-T 8:00 AM to 3:30 PM



## About This Kit

Superoxide dismutases (SODs) play a major role in antioxidant defense regulations by catalyzing the dismutation of superoxide anion free radical ( $O_2^{\cdot-}$ ) into molecular oxygen and hydrogen peroxide ( $H_2O_2$ ). Therefore, the amount of enzyme in cellular and extracellular environment is crucial for the prevention of oxidative stress damages.

TPR's SOD Assay Kit is a well-established tool for measuring SOD in **serum, plasma, tissue, cell lysates and other biological fluids**. This kit measures SOD activity by utilizing tetrazolium salt, which produces a water-soluble formazan dye upon reduction with superoxide anion (Fig. 1). The formazan formation rate is inhibited by presence of SOD in environments and is measurable photometrically.




Fig. 1: Scheme of the SOD assay

## Performance

- **Precision:** Human SOD with 8 replicates on 3 different days showed the intra and inter assay coefficient of variation 8.0% and 7.1% respectively.
- **Application:** This kit measures the activity of all isoforms of SOD enzyme and should work with all SODs (Cu/Zn, Mn and Fe SOD).

## Components

Item Label	Item	48 Tests Kit	96 Tests Kit
Reagent 1 (R1)	Chromogenic Reagent	1 ml	1 ml (2 vials)
Reagent 2 (R2)	SOD Enzyme Solution	15 $\mu$ l	30 $\mu$ l
Reagent 3 (R3)	SOD Assay Buffer	20 ml	20 ml (2 vials)
Reagent 4 (R4)	Sample Dilution Buffer	12 ml	12 ml
-	96 Well Microplate	1 Plate	2 Plates
-	Technical Manual	1 Manual	1 Manual

## Storage and stability

- This kit is stable for up to 6 months from date of receipt.
- Components are stable for 2 months after preparations, except Ready R2 Solution which is stable for 3 weeks only.

## Technical Service Contact Information

Phone: +98 21 86703137

Email: tpr.assist@gmail.com

Hours: S-T 8:00 AM to 3:30 PM



## Required Materials Not Provided

- Distilled water
- Disposables tubes.
- PBS or 150 mM KCl (if using tissue samples)
- Dounce homogenizer (if using tissue)

## Required Instrumentation

- Multiwell micro-plate reader (capable of measuring absorbance between 430-450 nm).
- Microfuge
- Incubator (37°C)
- 2-20  $\mu$ l and 20-200  $\mu$ l pipettes (multi-channel pipettes are preferred)

## Preparation of Sample Solutions

### ➤ Important Notes:

- ✓ All samples should be assayed immediately upon extraction or stored at -80°C for up to 1-2 months.
- ✓ Always thaw samples on ice.

### Cell Samples:

1. Harvest the amount of cells necessary for each assay ( $2 \times 10^6$  cells/ml recommended). For harvesting adherent cells, don't use proteolytic enzymes (like trypsin); rather use a rubber policeman.
2. Lyse cells by repeated cycles of freezing and thawing in ice cold PBS (pH 7.2-7.4).
3. Centrifuge at 14,000 x g for 5 minutes at +4°C.
4. Collect the supernatant and transfer to fresh tubes and place on ice.

### Tissue Samples:

1. Weight out the amount of tissue necessary for each assay (10-100 mg wet tissue).
2. Perfuse or rinse the tissue with PBS or KCl (150 mM), to remove any red blood cells and clots.
3. Homogenize the tissue in 1 ml ice cold PBS or 150 mM KCl.
4. Centrifuge at 14,000 x g for 5 minutes at +4°C.
5. Collect the supernatant and transfer to fresh tubes and place on ice.

**Note:** RIPA buffer (without SDS) supplemented with PMSF protease inhibitors can also be used as lysis buffer.

**Note:** If it is desired to measure SOD activity from cytosol and mitochondria separately, tissue/cell samples can be prepared according to method described by *Mattiazzi et al (2002)*. *JBC 277: 29626-33*.

### Plasma and Blood Samples:

1. Collect whole blood using an anticoagulant such as heparin, EDTA, or citrate.
2. Centrifuge at 1,000 x g for 10 minutes at +4°C.
3. Carefully pipette off the plasma layer without disturbing the buffy layer, transfer to fresh tubes and place on ice.

**Note:** Plasma can be diluted 3-10 times with Sample Buffer prior to assay SOD activity.

### Technical Service Contact Information

**Phone:** +98 21 86703137

**Email:** tpr.assist@gmail.com

**Hours:** S-T 8:00 AM to 3:30 PM



4. Remove the buffy layer from the red cell pellet.
5. Re-suspend the erythrocytes in four times its volume of ice cold distilled water.
6. Centrifuge at 10,000 x g for 15 minutes at +4°C.
7. Collect the supernatant (erythrocyte lysate) and transfer to appropriate size tubes and place on ice.

**Note:** Red Blood Cell lysate can be diluted 100 times with Sample Buffer prior to assay SOD activity.

### Serum Samples

1. Collect whole blood without using an anticoagulant.
2. Allow blood to clot for 30 minutes at 25°C.
3. Centrifuge the blood at 2000 x g for 15 minutes at 4°C and
4. Carefully pipette off the yellow serum layer without disturbing the buffy layer and transfer to new tubes and place on ice.

**Note:** Serum can be diluted 3-10 times with Sample Buffer prior to assay SOD activity.

### Reagent Preparation

- 1) **Working Solution:** Mix 1 mL of R1 with 19 mL of R3. Aliquot working solution so that you have enough volume to perform the desired number of assays. Working solution is stable for 2 months if stored at 4°C in dark.
- 2) **R2 Ready Reagent:** Spin R2 Reagent for 30 seconds and mixed well with pipette before dilution. Dilute 15 µL of R2 with 600 µL of R4 (48 Tests Kit) / Dilute 30 µL of R2 with 1.2 ml of R4 (96 Tests kit). R2 Ready Reagent is stable up to 3 weeks at 4°C.

### Things to Note

- For comparing SOD activities between different samples, these can be normalized on protein concentration or cell number. If necessary, the provided assay buffer can be used for diluting samples.
- All prepared reagents and samples must be equilibrated to room temperature just prior to use, and gently agitate.
- Assay all controls and samples in duplicate if possible.
- The SOD inhibition ratio can arrive 100%.
- Before the formal experiment, it needs to choose one or two samples for diluting a series of diluent and determine the dilution factor when the SOD inhibition ratio is 40%-60%.
- Since the superoxide will be released immediately after the addition of working solution, use a multiple channel pipette to avoid reaction time lag of each well.

## Technical Service Contact Information

Phone: +98 21 86703137

Email: tpr.assist@gmail.com

Hours: S-T 8:00 AM to 3:30 PM



## Assay Protocol

As described below, add samples and reagents in each well.

Blank 1 = 10  $\mu$ L ddH<sub>2</sub>O

Blank 2 = 10  $\mu$ L sample

Blank 3 = 10  $\mu$ L ddH<sub>2</sub>O

Sample wells = 10  $\mu$ L sample

Component	Sample ( $\mu$ L)	Blank 1 ( $\mu$ L)	Blank 2 ( $\mu$ L)	Blank 3 ( $\mu$ L)
Sample solution	10	0	10	0
ddH <sub>2</sub> O	0	10	0	10
Working solution	200	200	200	200
R2 ready reagent	10	10	0	0
R4	0	0	10	10

1. Add 200  $\mu$ L of Working Solution into each well.
2. Add 10  $\mu$ L of R4 to Blank 2 and Blank 3.
3. Add 10  $\mu$ L of R2 ready reagent to each sample well and Blank1.
4. Shake for 15 seconds and incubate at 37°C for 20 minutes.
5. Read the absorbance at 440-460 nm using a plate reader.

## Calculations

**Definition:** when SOD inhibition ratio in this reaction system reaches 50%, the corresponding enzyme level is 1 SOD activity unit (U).

### Formula for serum/plasma samples:

$$\text{Inhibition ratio of SOD (\%)} = \frac{(A_{\text{blank1}} - A_{\text{blank3}}) - (A_{\text{sample}} - A_{\text{blank2}})}{(A_{\text{blank1}} - A_{\text{blank3}})} \times 100$$

A = Absorbance

$$\text{SOD activity (U/mL)} = \text{Inhibition ratio of SOD (\%)} \div 50\% \times \left(\frac{220 \mu\text{l}}{10 \mu\text{l}}\right) \times \text{dilution factor of sample}$$

### Formula for tissue and cells:

$$\text{Inhibition ratio of SOD (\%)} = \frac{(A_{\text{blank1}} - A_{\text{blank3}}) - (A_{\text{sample}} - A_{\text{blank2}})}{(A_{\text{blank1}} - A_{\text{blank3}})} \times 100$$

A = Absorbance

$$\text{SOD activity (U/mg prot)} = \text{Inhibition ratio of SOD (\%)} \div 50\% \times \left(\frac{220 \mu\text{l}}{10 \mu\text{l}}\right) \times \text{dilution factor of sample} \div \text{Protein concentration of sample (mg prot/mL)}$$

### Technical Service Contact Information

**Phone:** +98 21 86703137

**Email:** tpr.assist@gmail.com

**Hours:** S-T 8:00 AM to 3:30 PM



### References

- Kussmaul L, Hirst J. 2006. The mechanism of superoxide production by NADH: ubiquinone oxidoreductase (complex I) from bovine heart mitochondria. *PNAS* **103**: 7607-12.
- Passos JF, Von Zglinicki T. 2006. Oxygen free radicals in cell senescence: are they signal transducers? *Free Radic Res* **40**:1277-83.
- Mattiazzi, M., D'Aurelio, M., Gajewski, CD., et al. 2002. Mutated human SOD1 causes dysfunction of oxidative phosphorylation in mitochondria of transgenic mice. *J biological chem.* **277(33)**:29626-33.
- Szeto HH. 2006. Mitochondria-targeted peptide antioxidants: novel neuroprotective agents. *AAPS J* **8**:E521-31.
- Robak J, Gryglewski RJ. 1988. Flavonoids are scavengers of superoxide anions. *Biochem Pharmacol* **37**:837-41.
- Nebot C, Moutet M, Huet P, Xu J-Z, Yadan J-C, Chaudiere J. 1993. Spectrophotometric assay of superoxide dismutase activity based on the activated auto-oxidation of a tetracyclic catechol. *Anal biochem* **214**:442-51.
- Okado-Matsumoto A, Fridovich I. 2001. Subcellular distribution of superoxide dismutases (SOD) in rat liver. *J Biol Chem* **276**:38388-93.